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
**Sensitive determination of elements in lithium batteries using the Thermo Scientific iCAP PRO XP ICP-OES**

*Authors: Jingliang He, Xiaobo Li, Fei Wang, Ming Jiang, and Jianping Cui*  
Thermo Fisher Scientific, Shanghai, China  
Thermo Fisher Scientific, Bremen, Germany

**Keywords:** Lithium battery, iCAP PRO XP ICP-OES, cathode material, battery materials, fast detection, robust analysis, quantification

**Goal**  
This note describes a fast analytical method for the determination of major and trace elements in the battery cathode material in lithium-ion batteries using the Thermo Scientific™ iCAP™ PRO XP ICP-OES.

**Introduction**  
The continuous development of lithium ion batteries and the research that relates to the behavior of the energy sector as it moves away from fossil fuels. To improve the health of production, the Chinese national standard GB 31137-2014 was established. All new lithium-ion power products must meet the requirements of GB 31137-2014. The most critical of these requirements include lithium, cobalt, nickel, and manganese, and potassium carbonate, as the cathode material in most cases. Lithium batteries only rarely release an cathode material that gradually releases active material through batteries. When used batteries and lithium ion phosphate batteries, there is a high volatility of graphite.

**Determination of elemental impurities in graphite powder for lithium-ion battery anodes**

*Author: Thermo Fisher Scientific Inc.*  
**Keywords:** Lithium ion batteries, elemental impurities, ICP-OES, graphite, anode

**Goal**  
To develop a robust and reliable method for the determination of elemental impurities in graphite powder samples for lithium-ion battery anodes, using the Thermo Scientific™ iCAP™ PRO XP ICP-OES. See introduction. This note will demonstrate the performance of the method with respect to detection limits (DL) of 10 ppm for major impurities (Ag, Pb, Fe, Ni, Cu) and quantification of these elements. It will also illustrate the recovery of the method developed through spike recovery analysis of the in-house spiked samples.

**Introduction**  
For lithium ion batteries, the negative electrode (anode) material is generally made from graphite powder. Graphite powder is made by the separation of graphite produced in an early method, chemically exfoliation, and non-metallic material with good mechanical properties and high-temperature oxidation, and corrosion resistance. It also has a large surface area and is often used as a high lithium ion reaction capacity and does not charge volume with reaction of structure. In addition, graphite powder can be recycled through reuse oxidation and purification processes to generate a core shell structure that can improve the energy-discharging performance and increase the anode lifetime. Graphite powder has become the main lithium ion battery anode material in use today in smaller consumer goods, such as mobile phones, as well as in electrical vehicles.

Through optimization of the sample pre-treatment process and dilution operation, together with rigorous selection of the most suitable anode anionizer, a robust, accurate and reliable ICP-OES analytical method for the determination of 13 trace elements in graphite powder materials for lithium ion batteries was developed.

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- IEC/TC 35 : Primary cells and batteries
- IEC/TC21 & TC69 : Secondary cells and batteries
- IEC/SC21A : Secondary cells and batteries containing alkaline or other non-acid electrolytes

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- ISO/WD 10655 : the determination of metal ions content of lithium hexafluorophosphate (LiPF<sub>6</sub>) by inductively coupled plasma optical emission spectrometry (ICP-OES).
- ISO/AWI 16398 : Lithium chloride - Determination of impurities by ICP-OES
- GB/T 26300-2020 : Nickel Cobalt Manganese composite hydroxide