

# MASSARRAY-BASED METHOD FOR DETECTION AND IDENTIFICATION OF FOODBORNE BACTERIA



**Dr. Namfon Suebwongsa**  
Product specialist  
Lifomics Co., Ltd.

# Contents

1

## MassARRAY Technology

- Introduction, MassARRAY principle
- System workflow
- Sample process journey
- Report & data analysis

2

## Foodborne Pathogen Detection

- Rationale and objective
- Sensitivity and specificity
- Limit of detection (LOD)
- Field sample detection

3


## Summary

4

## Q&A and Discussion

# KEY MESSAGE


เหนื่อยไหมกับการเลี้ยงเชื้อจำนวนมากเพื่อตรวจสอบการปนเปื้อนของจุลชีพก่อโรคในสิ่งส่งตรวจ??



การใช้เทคโนโลยีการตรวจ SNPs เพื่อการตรวจเชื้อก่อโรค



# PRINCIPLE OF MASSARRAY TECHNOLOGY




Cytosine (C)	Adenine (A)	Guanine (G)	Thymine (T)
247.2 Da	271.2 Da	287.2 Da	327.1 Da

Multiplex PCR





MALDI-TOF  
Mass Spectrometry

MassARRAY Detection Range  
4,000 – 9,000 Daltons






# SYSTEM WORKFLOW




# SAMPLE PROCESSES JOURNEY



## PCR/SAP/Extension


Multiplex end-point PCR followed by a single base extension reaction

**DNA SAMPLING 5-10 ng/ well**



## Target Sequence Detection

Automated analyte transfer and data acquisition



## Data Analysis

Data display and report generation

96 FORMAT



Run Time: 5.5 Hrs\*  
Hands-on Time: 25 Min



Run Time: 100 Mins\*  
Hands-on Time: 5 Min



Hands-on Time:  
5-10 Min



# MassARRAY: A High-Throughput Solution for Rapid Detection of Foodborne Pathogens in Real-World Settings

Namfon Suebwongsa<sup>1\*</sup>, Surasak Jiemsup<sup>1</sup>, Pannita Santiyanont<sup>1</sup>, Piyapha Hirunpatrawong<sup>2</sup>, Pornsiri Assawapirin<sup>2</sup>, Monthatip Thongkum<sup>2</sup>, Prakaymars Panumars<sup>2</sup>, Nipa Chokesajjawatee<sup>1</sup>, Supaporn Wongsrichai<sup>3</sup>, Pichet Koomba<sup>3</sup>, Suganya Yongkiettrakul<sup>1\*</sup>

<sup>1</sup>National Center for Genetic Engineering and Biotechnology (BIOTEC), Thailand, <sup>2</sup>Lifomics Co., Ltd., Bang Khae, Bangkok, THAILAND, Thailand, <sup>3</sup>Department of Livestock Development (Thailand), Thailand





# FOODBOURNE DISEASES



## Foodborne diseases in the WHO South-East Asia Region

Every year

 **>150 million** people fall ill } **>175 000** die

 including **60 million** children <5 fall ill } **50 000** die

Diarrhoeal diseases are responsible for majority of deaths. Key causes:

-  Norovirus
-  Non-typhoidal *Salmonella*
-  Pathogenic *E. coli*



Region has **>1/2** the people globally who are infected and die from **typhoid fever or hepatitis A**

**FOODBORNE DISEASES ARE PREVENTABLE. EVERYONE HAS A ROLE TO PLAY.**

Every year foodborne diseases cause:

almost **in 10** people to fall ill

**33 million** healthy life years lost

Foodborne diseases can be deadly, especially in children <5

 **420 000** deaths

 Children account for almost **1/3** of deaths from foodborne diseases

For more information: [www.who.int/foodsafety](http://www.who.int/foodsafety)

**#SafeFood**

Source: WHO Estimates of the Global Burden of Foodborne Diseases. 2015.



**World Health Organization**



# FOOD SUPPLY CHAIN



## Food supply chain




## Food safety monitoring and regulation



Food safety legislation, regulation and standards


Other legislation, regulation and standards




# OBJECTIVE

To develop an efficient **high-throughput assay** for bacterial detection and identification based on **MassARRAY<sup>®</sup> technology**


- *Reliable and feasible*
- *Non labor intensive*
- *Less time consuming*
- *More cost effectiveness*





**BIOTEC**  
a member of NSTDA


# BACTERIAL FOODBORNE PANEL



Campylobacter Clostridium perfringens Listeria monocytogenes E. coli Staphylococcus aureus Salmonella Shigella

## 10 BACTERIAL FOODBORNE PATHOGENS IN A SINGLE REACTION

- Salmonella* spp.
- Enterococcus faecalis*
- Campylobacter coli*
- Enterococcus faecium*
- Campylobacter jejuni*
- Staphylococcus aureus*
- Listeria monocytogenes*
- Clostridium perfringens*
- Escherichia coli* and *Shigella* spp.




$$*\% \text{Sensitivity} = \frac{\text{true positive}}{\text{true positive} + \text{false negative}} \times 100$$

$$**\% \text{Specificity} = \frac{\text{true negative}}{\text{true negative} + \text{false negative}} \times 100$$



# MassARRAY-based bacterial detection assay

Detection of 10 bacterial foodborne pathogens in a single reaction




*Enterococcus faecalis*  
*Enterococcus faecium*



*Clostridium perfringens*




*Staphylococcus aureus*




*Listeria monocytogenes*




*Escherichia coli*



*Shigella* spp.



*Salmonella enterica*



*Campylobacter jejuni*  
*Campylobacter coli*

# Research workflow

## 1. Target gene and Primer design

- F/R primer
- Extension primer

## 2. Assay specificity testing

- Reference gDNA
- Laboratory isolates

## 3. Assay validation

- 103 field samples (meat)



# Primers and Assay designed

**Table 1.** Primers used in this study

Targeted species	Primers	PCR primers (forward, reverse)		PCR amplicon size (bp)	Extension primers			Mass of SBE products (Da)
		Tm (°C)	% GC		Tm (°C)	% GC	Mass (Da)	
Bacteria	Bac16 1-1	62.8, 66.9	54.5, 48.4	105	45.8	38.9	5545.6	5816.8
<i>Campylobacter coli</i>	Camp005	63.5, 61.9	50.0, 45.2	150	50.0	42.9	6747.4	7034.6
	Camp006	63.6, 59.8	50.0, 43.3	105	49.6	39.1	7511.9	7783.1
<i>Campylobacter jejuni</i>	Camp002	64.1, 63.5	50.0, 50.0	108	46.5	27.3	6781.5	7108.5
<i>Clostridium perfringens</i>	Clos001	61.7, 63.8	46.7, 50.0	101	46.2	21.7	7078.7	7365.9
<i>Escherichia coli/Shigella</i> spp.	Eco001N	62.0, 63.9	46.7, 50.0	108	48.6	36.8	5795.8	6043.0
<i>Enterococcus faecium</i>	Ent001	61.2, 63.9	46.7, 50.0	99	48.4	28.0	7670.0	7957.2
<i>Enterococcus faecalis</i>	Ent003	62.6, 63.9	46.7, 50.0	141	49.3	50.0	5425.5	5752.6
<i>Listeria monocytogenes</i>	Lis001	62.5, 61.8	50.0, 46.7	99	45.7	56.2	4955.2	5282.3
	LisG	-	-	99	52.1	62.5	4979.2	5306.3
<i>Salmonella</i> spp.	Sal002	63.4, 63.8	50.0, 50.0	106	51.1	42.9	6365.2	6652.4
<i>Staphylococcus aureus</i>	Stap001	58.6, 59.3	38.7, 43.3	94	46.0	35.0	6238.1	6565.2

**Table 2.** Genomic DNA samples used as reference DNA templates in this study.

	No.	Samples (gDNA)	Strains	Sources
Bacteria	1	<i>Campylobacter coli</i>	JV20	ATCC
	2	<i>C. jejuni</i>	D3071	ATCC
	3	<i>Clostridium perfringens</i>	WAL-14572	ATCC
	4	<i>Escherichia coli</i>	B171	ATCC
	5	<i>Enterococcus faecalis</i>	TX0104	ATCC
	6	<i>E. faecalis</i>	TUSoD Ef11	ATCC
	7	<i>Enterococcus faecium</i>	TX0133a04	ATCC
	8	<i>Listeria marthii</i>	FSL S4-120	ATCC
	9	<i>L. monocytogenes</i>	F6900	ATCC
	10	<i>L. monocytogenes</i>	Li21	ATCC
	11	<i>Salmonella enterica</i>	15/5	ATCC
	12	<i>S. enterica</i>	LT2	ATCC
	13	<i>Shigella</i> spp.	NMICID504	Human <sup>a</sup>
	14	<i>Staphylococcus aureus</i>	CM05	ATCC
Protozoa	15	<i>Crithidia fasciculata</i>	ND	Human <sup>b</sup>
	16	<i>Leishmania martiniquensis</i>	ND	Human <sup>b</sup>
Virus	17	Severe acute respiratory syndrome coronavirus 2	ND	Human <sup>b</sup>
Human	18	Human	No.4312660	Applied Biosystems, California, United States.
Plant	19	Plant	ND	Teak leaf <sup>c</sup>



# Validation and evaluation of MassARRAY-based assay

## Determination of sensitivity and specificity of MassARRAY-based assay with gDNA

**Table 3.** Bacterial isolates used

Samples (gDNA)	Reference (no.)	Lab isolate (no.)	Total (no.)
<i>Staphylococcus aureus</i>	2	7	9
<i>Salmonella spp.</i>	5	13	18
<i>Shigella spp./Escherichia coli</i>	2	6	8
<i>Listeria monocytogenes</i>	4	4	8
<i>Enterococcus faecium</i>	1	4	5
<i>Enterococcus faecalis</i>	1	10	11
<i>Campylobacter jejuni</i>	0	8	8
<i>Campylobacter coli</i>	0	3	3
<i>Clostridium perfringens</i>	0	5	5
Others species	3	7	10
Mixed templates	4	0	4
Total	22	67	89

# MassARRAY® System workflow



1

PCR AMPLIFICATION

2

SAP

3


SINGLE-BASE EXTENSION (SBE)

4

SAMPLE TRANSFER & DATA ECQUISITION

5

DATA ANALYSIS & REPORTING



15 MINS Hands-on Time

3 MINS Hands-on Time

5 MINS Hands-on Time

5 MINS Hands-on Time

# MassARRAY Data Generation and Interpretation

TypoAnalyzer SEQUENOM

File View Options Tools Help

Project Explorer Traffic Light Details

Customer/Project/Plat

	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B													
C													
D													
E													
F													
G													
H													

Chip

F0110992\_RECALL\_185

Chip Summary

Chip F0110992\_RECALL\_1856129\_ (1) (1)

Conservative Calls 170

Moderate Calls 2

Aggressive Calls 1

User Calls 0

Calls 173

No Calls 883

Total Possible Calls 1056

Assay

Bac16

Clos001

Cmp002

Cmp005

Cmp006

Ent001

Ent003

Lis001

Sal002

Stp001

Well: A01 Sample: Prim1-2\_1

Lis001

Intensity


Mass

Assay	Call	Description
Bac16	A	A.Conservative
Clos001		N.No-Alleles
Cmp002		N.No-Alleles
Cmp005		N.No-Alleles
Cmp006		N.No-Alleles
Ent001		N.No-Alleles
Ent003		N.No-Alleles
lacY_Ecoli08		N.No-Alleles
Lis001		D.Low Probability
Sal002		N.No-Alleles
Stp001	A	A.Conservative


Assay Info	Call Info	Peak Info	Area Info
Expected Peaks	Name	Sequence	Mass
Probe	UEP.Lis001	CAGGGAGAACATCTGG	4955.20
Analyte	T	CAGGGAGAACATCTGGA	5226.40
Analyte	A	CAGGGAGAACATCTGGT	5282.30

Automated Data ... Assay Summary Chip Summary Post Processing Clusters Call Cluster Plot Histogram Details






**Figure 1.** The chromatograms illustrate the specific molecular masses of SBE products corresponding to bacterial targets.




- (A) a negative control (DW).
- (B) SBE product of bacteria (Bac16 1–1).
- (C) SBE products of *C. jejuni*
- (D,E) SBE products of *C. coli*.
- (F) SBE products of *Cl. perfringens*
- (G) SBE products of *E. coli*
- (H) SBE product of *Shigella spp.*
- (I) SBE products of *E. faecium*.
- (J) SBE products of *E. faecalis*.
- (K) SBE products of *L. monocytogenes* variant 1.
- (L) SBE products of *L. monocytogenes* variant 2
- (M) SBE products of *Salmonella spp.*
- (N) SBE products of *S. aureus*.




**Figure 2.** The chromatograms representing negative results clearly illustrate specific molecular mass of the extension primers when human gDNA (A) and plant gDNA (B) were used templates in the MassARRAY-based assay. Similar chromatograms were observed for additional negative controls, including *Listeria marthii*, *Crithidia fasciculata*, *Leishmania martiniquensis*, and Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2). These chromatograms indicated consistent negative results across different non-targeted organisms, thus confirming the absence of cross-reaction in the MassARRAY-based assay system.

### *Salmonella* spp./*C. jejuni*



- ➔ 16s rRNA
- ➔ *Salmonella* spp.
- ➔ *C. jejuni*
- ➔ *E. faecium*

### *Salmonella* spp./*C. jejuni*/*E. faecium*



**Figure 3.** The chromatograms representing negative results clearly illustrate specific molecular mass of the extension primers when using mixture template gDNA

# Sensitivity and Specificity of MassARRAY-based assay

**Assay sensitivity 98.8%**

**Assay specificity 100 %**

Sample (genomic DNA)	Number of tested samples		Total	Base calling			%Sensitivity	%Specificity
	Reference strains	Lab isolated strains		Detected	Not detected	Cross-reaction		
<i>Staphylococcus aureus</i>	2	7	9	9	0	0	100	100
<i>Salmonella</i> spp.	5	13	18	18	0	0	100	100
<i>Shigella</i> spp. / <i>E. coli</i>	2	6	8	8	0	0	100	100
<i>Listeria monocytogenes</i>	4	4	8	8	0	0	100	100
<i>Enterococcus faecium</i>	1	4	5	5	0	0	100	100
<i>Enterococcus faecalis</i>	1	10	11	11	0	0	100	100
<i>Campylobacter jejuni</i>	0	8	8	8	0	0	100	100
<i>Campylobacter coli</i>	0	3	3	3	0	0	100	100
<i>Clostridium perfringens</i>	0	5	5	5	0	0	100	100
Others species								
<i>Clostridium</i> spp.	0	4	4	3 <sup>a</sup>	1	0	75	100
<i>Klebsiella pneumoniae</i>	0	3	3	3 <sup>a</sup>	0	0	100	100
<i>Listeria innocua</i>	2	0	2	2 <sup>a</sup>	0	0	100	100
<i>Campylobacter lari</i>	1	0	1	1 <sup>a</sup>	0	0	100	100
Total	18	67	85	84	1	0	98.8	100

Data summarized from two independent experiments.

<sup>a</sup> Base calling for *16s rRNA* gene only.


# Validation and evaluation of MassARRAY-based assay

## Determination of limit of detection (LOD) with bacterial cells

Bacterial culture incubated at 37° C or 41° C for 18 h



Adjusted OD 600 = 1.85 and collected bacterial cells



100 µL for gDNA extraction using QuickExtract™




MassARRAY detection



# LOD of MassARRAY-based assay

- LOD  $\sim 10^2 - 10^3$  cells of *E. coli*, *C. perfringens*, *C. jejuni*, and *Salmonella* spp.
- LOD  $\sim 10^4 - 10^5$  cells of *E. faecalis*, *S. aureus*, *C. coli*, *L. monocytogenes*, and *E. faecium*




Data obtained from three independent experiments.

# Validation of MassARRAY with unknown specimens

In collaboration with  
Bureau of Quality Control of Livestock  
Products  
(สำนักงานตรวจสอบคุณภาพสินค้าปศุสัตว์ (สตส.) กรมปศุสัตว์)

# Pathogen detection procedure of Bureau of Quality Control of Livestock Products



# Lab Workflow

Collected 5-10 mL, centrifuged at 4,000-7,000 rpm 10 min.  
and decarded media.



The DNA extraction procedure employed the  
MagPurix® Bacterial DNA Extraction Kit (ZP02006),




1 µl of gDNA for PCR by using EconoTaq® Polymerase



MassARRAY



# Pathogen detection procedure of Bureau of Quality Control of Livestock Products



# Result: Validation of MassARRAY with unknown specimens

Sample no.	Sample ID	Bacterial culture-based method					MassARRAY-based assay				
		Bacteria	<i>E. coli/Shigella</i> spp.	<i>Enterococcus</i> spp.	<i>Salmonella</i> spp.	<i>S. aureus</i>	Bacteria	<i>E. coli/Shigella</i> spp.	<i>Enterococcus</i> spp.	<i>Salmonella</i> spp.	<i>S. aureus</i>
1	4,627	P	N	N	N	N	P	N	N	N	N
2	4,628	P	N	N	N	N	P	N	N	N	N
3	4,646	P	P	N	N	N	P	P	N	N	N
4	5,543	P	P	N	N	P	P	P	N	N	P
5	5,545	P	N	N	N	N	P	N	N	N	N
6	5,546	P	N	N	N	N	P	N	N	N	N
7	5,555	P	N	N	N	N	P	N	N	N	N
8	5,556	P	P	P	N	N	P	P	P*	N	N
9	5,571	P	P	N	N	N	P	P	N	N	N
10	38,532	P	P	P	P	N	P	P	P*	P	N
11	38,633	P	P	P	P	P	P	P	P*	P	P
12	38,809	P	P	P	N	N	P	P	P*	N	N

N, negative result; P, positive result. \*The method specifically identified the species as *E. faecalis*.

91 samples showed positive bacteria  
 12 samples presented positive target bacteria  
 100% concordance with standard culture method

# SUMMARY

- The developed MassARRAY-based assay demonstrated its high-throughput performance in simultaneously identifying targeted species in a single assay with high specificity.
- Its successful application to field samples positioned it as a compelling choice for the rapid detection of foodborne bacterial pathogens in real-world field settings.
- Furthermore, the cost effectiveness, time efficiency, and accuracy inherent to MassARRAY technology make it highly suitable for adoption in food industrial laboratories, particularly those engaged in quality control of food products.

# MassARRAY-based application for food safety and security



# Acknowledgement



**National Center for Genetic Engineering and  
Biotechnology**  
National Science and Technology Development Agency



**Bureau of Quality Control of Livestock Products, Department of Livestock  
Development, Ministry of Agriculture and Cooperatives, Pathum Thani,  
Thailand**





**BIOTEC**<sup>1</sup>  
a member of **NSTDA**



---

**Thank you**

---

**Q&A**  
**Discussion**

# Efficiency of MassARRAY detection

Sample	Sample type	Number
Meat sample	DF/BPB	103 (day 0)
Total 103 samples	BPW	103 (day 1)
	RVS /MKTTn broth	103 (Day 2)
	Selective plates	11 (Day 3)*

\* The day of sample obtained from Bureau of Quality Control of Livestock Products.